



# Z-Dis

## Dissolution Monitor

Release 2.09  
for Windows NT 4.0

## Software - Manual

English

---

Document created: 18.12.96 by J. Mondry  
Last update: 10.07.98 by J. Mondry

© 96-98 by Carl Zeiss Jena GmbH, OEM/M2 Process Sensorics Dept.



〒164-0011  
東京都中野区中央4-4-5第一小林ビル  
Tel: 03-5328-2858 Fax: 03-5328-2859

# Contents

**CONTENTS..... 2**

**1. HARDWARE AND SOFTWARE SYSTEM SPECIFICATION..... 3**

**2. INSTALLATION AND INI FILES ..... 4**

    INSTALLATION ..... 4

    DISSTEST.INI..... 5

    ZY.INI..... 8

**3. DATA STORAGE AND FILE SYSTEM..... 9**

**4. THE AUTOMATIC MODULE..... 11**

**5. THE MANUAL MODULE..... 15**

    Z-DIS MUX PARAMETER FILE FORMAT ..... 18

**6. RELEASE HISTORY ..... 19**

**SUPPLEMENTS..... 24**

    ALGORITHM FOR SMOOTHING AND DERIVATIVES ACCORDING TO SAVITZKY-GOLAY ..... 25

    LABVIEW GRAPH HANDLES ..... 26

    ZEISS ASPECT PLUS DATA FILE FORMAT ..... 28

---

# 1. Hardware and software system specification

The Z-Dis Dissolution Monitor Software is designed and compiled using exclusively National Instruments™ LabView™ development system (current release 4.01). The software is programmed under and designed for MS Windows NT 4.0. The software is not portable to other operating systems.

Being a Win NT 4.0 application, the Z-Dis Dissolution Monitor Software requires the following minimum hardware and software environment:

- IBM compatible Pentium™ 133 with at least 32 MB RAM (64 MB recommended)
- about 50 MB of free harddisk capacity
- SVGA graphics card
- Mouse
  
- MS Windows NT 4.0 (logged on to a NT network)
- NT graphics mode set to 1024 x 768 x 256 colors, small fonts (!)
- A suitable printer (Laserjet type) and a suitable storage device (harddisk) must be logged on to the network for remote protocol printout and data storage. The machine running Z-Dis must have access to the network printer and storage device.
  
- No Eco-PC (Power savings etc.) settings installed or activated (\*)  
*Because **Z-DIS** is designed for long duration runs it is strongly recommended to uninstall or disable, respectively, any Eco-PC settings of your PC ( the likes of e.g. processor shut-down if no input is made by mouse or keyboard etc.)*
  
- The ZYMARK automatic system must be properly connected to the I/O ports at the back of the ZEISS MCS spectrometer.

*For a pinout of the MCS 500 I/O ports see Supplement A.*

---

## 2. Installation and INI files

### *Installation*

This release is delivered complete with a standard installation routine.

To start the installation:

- insert disk 1 into the floppy drive of your PC
- run ZDISnnn.EXE from the Windows file manager

The software will now be installed, the required directories etc. will be automatically created. The default installation directory is **Z-Dis\_Robotic\_NT**.

A program manager group will be created holding several program icons.

The software package includes two executables:

- **Disstest.EXE** is the Automatic Module to perform the dissolution runs and communicate with the ZYMARK robotic system. To start the automatic module double-click on the AUTOMATIC icon in the program manager group.
- **Manual.EXE** is the Manual Module to create parameter sets for the Automatic Module. To start the manual module double-click on the MANUAL icon in the program manager group.

You can also view and edit DISSTEST.INI (Z-DIS Initialisation file), DEFAULT.PAR (a set of default measurement parameters) and VERSION.TXT (info & release history) from the program manager group.

## ***DISSTEST.INI***

The initialisation file DISSTEST.INI holds essential information for the program execution. It has to be in the applications main directory. Both the automatic as well as the manual modules are referring to this INI file.

Example INI file:

```
[General]
Version=Version 2.01 NT
Language=1
Debug=1

[Files]
CurDataFile=\\Jeeves\Pfizer\Zymark\Data\Test.zrd
BackupDataFile=d:\Development\Projects\Z-Dis_Robotic_NT\TempReportFile.zrd
CurParaFile=d:\Development\Projects\Z-Dis_Robotic_NT\Parameters\default.PAR

[MCS500]
Simulation=yes
MCSType=MCS 551 UV
NrOfDevices=1
MUX=No

[Dissolution]
SimMode=2
RobotINI=d:\Development\Projects\Z-Dis_Robotic_NT\ZY.INI
StoreRawData=0

[DISSTEST]
FPFont="0" 16
BDFont="0" 16
postScriptLevel2=False
colorPrinting=False
IconFont="0" 16
simpleDiagramHelp=False
appFont="0" 13
useDefaultTimer=False
noBoxAroundName=True
```

Some of the parameters in DISSTEST.INI can be manually modified by the experienced user. In the following description these are marked by an (\*).

### **[General]**

Version=Version 2.01 NT

Language=1 English

Debug=0 internal

## [Files]

CurDataFile=	last used ZYMARK data log file
e.g.: = \\Jeeves\Pfizer\Zyemark\Data\Test.zrd	
BackupDataFile=	Backup file for data storage in case of network failure
e.g.: d:\Development\Projects\Z-Dis_Robotic_NT\TempReportFile.zrd	
CurParaFile=	last used parameter file
e.g.: d:\Development\Projects\Z-Dis_Robotic_NT\Parameters\default.PAR	

## [MCS500]

Simulation=no	<b>Spectrometer simulation mode (*)</b>
no	no simulation, actual measurement
results	spectrometer on-line: -> all port commands and measurements are carried out, but spectra will not be evaluated for concentration, color etc.; instead dummy values will be generated
yes	no spectrometer(s) on/line, all spectra and eval's will be emulated
MCSType=MCS 551 UV	Type of Zeiss spectrometer used
NrOfDevices=1	Number of spectrometers in the pipeline
MUX=No	Shutter or Multiplexer device ?

## [Dissolution]

SimMode=0	<b>ZYMARK system simulation mode (*)</b>
0	robot on-line
1	use manual DIO-Emulation
2	automatic simulation mode
RobotINI=	path and name of the ZYMARK INI file
e.g.: = d:\Development\Projects\Z-Dis_Robotic_NT\ZY.INI	used for communication between ZEISS and ZYMARK programs.
StoreRawData=0	0/1 (*) specifies whether to store the raw spectral data of samples, background and standards. The files will be stored in ZEISS

---

ASPECTplus data file format and can be re-evaluated any time using the ASPECTplus software.

**[DISSTEST]**

This section hold various settings for the correct screen display, fonts etc. Don't change these items.

*For a description of the ZEISS ASPECT plus data file format see Supplement D.*

The INI file parameter [Dissolution] SimMode allows you to make use of an internal DIO simulation to manually trigger the single commands of a dissolution run. You will be prompted to confirm the next awaited command. Click on OK to confirm and continue.

## Zy.INI

The communication between the Zymark Software System and the Zeiss Z-Dis software is performed using a common file, Zy.INI, into which essential information for the next run is written by the Zymark system. These data will be read by Z-Dis upon the “re-read Zy.ini” command.

### Example of Zy.INI

```
[Zymark]
ZeissPARFile=d:\Development\Projects\Z-Dis Robotic NT\Parameters\default.PAR
ARDReference=123456-sd45
BatchSize=6
StartingVessel=1
VesselVol=500
SamplePoints=3
DosageSize=500
StdConcentration1=1.0
StdConcentration2=1.0
StdContainer1=1
StdContainer2=2
MaxBackgroundAbs=0.1
SystemSuitability=1
Ratio=5
Ratio2=5
RSD=5
QPoint=80
QPointSP=3
DataFilePresent=1
DataFilePath=\\Jeeves\Pfizer\Zymark\Data\test.zrd
```

<b>ZeissPARFile</b>	(Net-) Location of the parameter file to be used by Z-Dis
<b>ARDReference</b>	ARD Ref number to be included in result output strings
<b>BatchSize</b>	6 vessels
<b>StartingVessel</b>	1 or 6
<b>VesselVol</b>	[ml]
<b>SamplePoints</b>	Number of sample points
<b>DosageSize</b>	[mg/ml]
<b>StdConcentration1, StdConcentration2</b>	Concentration of the standards to be measured
<b>StdContainer1, StdContainer2</b>	Container for standard solution 1 and 2
<b>MaxBackgroundAbs</b>	Maximum tolerated background absorbance
<b>SystemSuitability</b>	Perform system suitability run before the run ? (0/1)
<b>Ratio, Ratio2, RSD</b>	Tolerances for Std1/2 ratio, Std 2/1 ratio and SST standard measurement RSD
<b>Qpoint, QPointSP</b>	Qpoint, and Qpoint-sample point for graphical display
<b>DataFilePresent</b>	Data file ready ? (0/1)
<b>DataFilePath</b>	(Net-) Location of Data file



### 3. Data Storage and File system

The Results generated during a run are stored in the common data file specified in the ZY.INI. The structure and contents of the individual lines provided by Z-Dis and the Zymark system are specified in detail in the Zymark documentation.

See Document: **PFZ 09 Dissolution System File Interface**

Author Paul Savage, Zymark

Excerpt from example data file:

```
...
PARAMETER_FILE_ZEISS,d:\Development\Projects\Z-Dis Robotic NT\Parameters\default.PAR,400.00,NORMAL,5,30,1
BLANK,0,Fri Aug 15 14:06:03 1997,0
...
...
STD2,1,2,Fri Aug 15 14:06:03 1997,1.0000,0.2904
STD2,2,2,Fri Aug 15 14:06:03 1997,1.0000,0.2906
STD2,3,2,Fri Aug 15 14:06:03 1997,1.0000,0.2897
STD2,4,2,Fri Aug 15 14:06:03 1997,1.0000,0.2897
STD2,5,2,Fri Aug 15 14:06:04 1997,1.0000,0.2891
STD2,6,2,Fri Aug 15 14:06:04 1997,1.0000,0.2891
STD2,7,2,Fri Aug 15 14:06:04 1997,1.0000,0.2892
STD1,1,1,Fri Aug 15 14:06:04 1997,1.0000,0.2899
SYSTEM_SUITABILITY,100.18,PASS,0.29,0.001,PASS,123456-sd45
...
...
BACKGROUND,1,Fri Aug 15 14:06:05 1997,0.00828
BACKGROUND,2,Fri Aug 15 14:06:05 1997,0.04950
BACKGROUND,3,Fri Aug 15 14:06:05 1997,0.08753
BACKGROUND,4,Fri Aug 15 14:06:05 1997,0.12224
BACKGROUND,5,Fri Aug 15 14:06:05 1997,0.15470
BACKGROUND,6,Fri Aug 15 14:06:06 1997,0.18468
...
```

If the Option StoreRawData is activated (Disstest.ini) all acquired spectra, i.e. background, blank, standard and sample readings, are stored into the data directory specified in the ZY.INI. The spectra are stored in ASPECT plus native spectra file format.

*For a description of the ZEISS ASPECT plus data file format see Supplement D.*

The file names are generated automatically by the Z-Dis program and indicate the contents of the file according to the following table.


### **Stored raw data files**

<b>Filename</b>	<b>Contents</b>
Yyyymmddrr Sample Vessel v.dat	[Samplepoints] Sample spectra
Yyyymmddrr Manual Mode Spectra.dat	Spectra generated in the manual mode
Yyyymmddrr Standard Spectra.dat	2 (or 8 in case of SST) Standard spectra
Yyyymmddrr Background Spectra.dat	6 Background spectra

*[yy..year, mm .. month, dd .. day, vv .. vessel nr., rr .. run nr ]*

## 4. The Automatic Module

Since the Z-Dis software runs in tight cooperation with the Zymark software system, it occupies only one half of the operator's screen.



---

Short description of screen elements.

**Dissolution Graph** displays the dissolution curves for any single or all of the dissolution vessels. The vessels to display can be specified in the checkbox array left of the dissolution graph.

**Time** shows the current duration of the dissolution run.

**Run** indicates the current run or method, resp.

**Cycle** indicates the current cycle or sample point, resp.

**Vessel** indicates the vessel that is currently being measured

**DIO In / Out** show the current I/O port status (communication with ZYMARK robot)

**Spectra Graph** Displays the current spectra obtained from the activated vessels. The upper selection box specifies whether **Background**, **Standard** or **Sample** spectra are displayed. The lower selection box specifies whether the spectra are displayed in **Energy**, **Absorbance** or **2<sup>nd</sup> deriv of Absorbance**.

All displays are continuously updated while the dissolution run is performed.

The Z-Dis Software is started and closed down by the Zymark software system.  
(A hidden function enables the F11 key to close down the Z-Dis Software)

The Automatic Module software awaits a well defined chain of commands from the ZYMARK robotic system. The structure of the command sequence is described in the following table.

ZYMARK command (DIO In)	ZEISS response (DIO Out)	Automatic Module software
		Reads initialisation parameters from ZDIS_MUX.INI, especially location of ZY.INI file and SIM mode
1	1, 4	Reads actualized ZY.INI, gets batch size, std conc. a.o. values, as well as ZEISS para file to use in this run; reads specified ZEISS para file, gets calibration wavelength etc.
2	1,2,4	Reads blank spectrum
3	1,3,4	System Suitability Test (optional): Read Std 1, compute and evaluate ratio, read 8 times Std 1 and compute and evaluate RSD.
3	1,3,4	Reads Std 2, computes preliminary calibration
2	1,2,4	Reads nn background spectra (nn.. nr. of vessels), checks absorbance at specified wavelength and stores background as reference for sample measurement. If one of the background absorbances exceeds MaxAllowedBackground the complete background cycle is to be repeated.
4	1,4	Starts the analysis timer
3	1,3,4	Reads nn * cc sample spectra (nr. of vessels * nr. of sample points). Each sample spectrum is corrected with the respective vessel's background reading. The concentration is determined from the preliminary calibration and displayed in the dissolution graph.
2	1,2,4	Reads second blank spectrum
3	1,3,4	Read first standard spectrum, computes mean calibration.

---


Now all the concentrations of the sample run are re-evaluated using the mean calibration. If RawDataStorage is ON, all collected spectra are stored in corresponding raw data files.

The whole sequence can now be repeated with the next method.

The results of the dissolution run will be added to the (remote drive) ZYMARK data log file during the run.

A more detailed description of the handshake communication between Zymark and Zeiss programs can be found in the Zymark system documentation.

## 5. The Manual Module



The Manual Module is used to create and evaluate parameter files for the automated dissolution runs. It allows you to observe the complete spectra and their 2nd derivatives in order to find the best calibration wavelength, integration time and other settings for the special analytical problem.

The filename and path of the current parameter file are displayed in the Headline Bar of this window.

*For a description of the graph manipulation handles (upper right of the graphs) see Supplement C.*

<b><i>Integration time</i></b>	Specifies the integration time of the MCS 500 diode array spectrometer.
<b><i>Average</i></b>	Specifies the number of spectra to be averaged by the spectrometer in each measurement.
<b><i>Standard Conc.</i></b>	Specifies the concentration of the standard that will be measured after activating the STANDARD button. The value is required to compute the concentration of a sample measurement.
<b><i>Calibration wavelength</i></b>	Specifies the wavelength at which to compute the calibration curve. The current calibration wavelength is indicated by a red line cursor in the graphs.
<b><i>Smoothing</i></b>	Determines the number of measuring points included in the smoothing algorithm for the computation of the 2nd derivative.
<b><i>Source</i></b>	Specifies whether the calibration curve is taken from the standard spectrum or from its 2nd derivative
<b><i>Sample Abs, Std Abs and % Dissolution</i></b>	Results of a sample measurement as computed from the current settings and blank/std readings.
<b><i>F1 Blank</i></b>	Blank measurement. The obtained spectrum will be stored as blank (reference) for standard and sample measurements.
<b><i>F2 Standard</i></b>	Standard measurement. The standard spectrum and its derivative will be displayed in yellow.
<b><i>F3 Sample</i></b>	Sample measurement. The obtained spectrum is blank-(reference-) corrected. The sample spectrum will be displayed. If a valid standard has been measured, the concentration is computed and displayed.
<b><i>F4 Load Para</i></b>	Allows the user to load a complete parameter set.



---

**F5 Save Para**

Allows the user to store the complete parameter set, i.e. calibration settings and spectrometer settings.

**F6 Load Spectra**

Allows the user to load a set of spectra.

**F7 Save Spectra**

Allows the user to store the current set of spectra.

**F8 Plot**

Plot current spectra and results.

**F10 Quit**

Quit Manual Module.

**Spectra Graph 1 & 2**

Displays measured blank, standard and sample spectra. The selection boxes specify whether the spectra are displayed in **Energy, Absorbance** or **2nd deriv** of Absorbance.

## ***Z-Dis MUX parameter file format***

A parameter file contains all necessary parameters for a complete dissolution run, i.e.:  
The parameter files are stored in plain ASCII (use notepad or any other ASCII viewer)

Example:

```
[GENERAL]
FileID=Z-Dis MUX Dissolution Monitor Parameter File
Time=21.11.1995/12:26:49
Origin=Dissolution Monitor V 1.54
Note=

[Spectrometer]
IntegrationTime=25
Accumulation=20

[Evaluation]
Smooth=8
CalibWavelength= 400.00
UseDeriv=0
```

## 6. Release History

### *Release 2.09 NT*

released : 04.08.98  
at Sandwich site

- changes according to Zymark FAX 15.7.98 and comments by S.Lobley of Pfizer
- std1, std2 & blank container read from ZY.INI correctly
  - std 2 8th reading now correctly numbered 8 in ZRD file
  - spectra file name now includes date & time to avoid overwriting in case of multiple system start at the same day
  - graph colors changed according to Win3.11 release
  - spectra stored with date time stamp in correct entry according to ASPECT file format spec
  - manual mode: graph plot now has correct scaling
  - std & sample reading indicators and conc indicator cleared upon blank or std measurement, resp.

All changes carried out according to Software Modification Procedure and according Software Modification reports have been generated.

### *Release 2.08 NT*

released : 24.02.98  
at Sandwich site

- started from 2.04 source code
- fixes according to 2.05 redone
- SST string carried over to subsequent runs
- new calibration cluster  
(std 2 data separated into SST(single & 6x) and final meas.)

### *Release 2.07 NT*

(rejected)  
released : 10.02.98

- BF in IO and Shutter bit logic
- updated firmware included (mcs501.btl of 1.9.97)
- new INI entry IOFlag now allows logic swap independently from shutter / MUX

```
example : [MCS500]      Simulation=No  
                      MCSType=MCS 551 UV  
                      NrOfDevices=1  
                      Dark=Shutter  
                      IOFlag=0
```

### *Release 2.06 NT*

(rejected)  
released : 04.02.98

changes according to FAX of 3.2.98 completed at Zeiss,  
replacement EXE file emailed to Zymark

- BF in SST data transfer to following runs

**Release 2.05 NT**

released : 28.11.97

changes according to FAX of 28.11.97 completed at Zeiss,  
only sent EXE file

- added absolute value in RSD evaluation (BF in RSD pass/fail check)
- changed ZY.ini entries to latest Zymark state, (BF in STD pass/fail check)

**Release 2.04 NT**

released at : 3.11.97

changes according to FAX of 29.10.97 completed at Pfizer Sandwich

- added [Directories] entry to INI file (spectra storage)
- checked output lines according to Zymark specification Release 1.08
- BF in vessel nr display at begin of run
- smaller fonts in Manual plot

**Release 2.03 NT**

released at : 6.10.97

changes according to FAX of 25.09.97

Automatic:

- All readings in std, backgr and sample report lines are now 7 decimal places
- System suitability data carried over from last SST, SST and Std calculation now correct
- clear screen immediately after final reading, then wait for switch 1
- run number updated upon re-reading zy.ini

Manual:

- BF in save & load spectra
- volume and dosage inputs removed, thus result is shown as conc.
- all results now 7 dp, input 7dp as well
- parameter default directory read from disstest ini

General:

- updated firmware file MCS501.BTL in use

**Release 2.02**

Released 02.09.97

All readings in std, backgr and sample report lines are now 6 decimal places  
QPoint line now is displayed as Q+5  
Dissolution graphs start at (0,0), times rounded to half minute

**Release 2.01**

Released 22.08.97

BF in SST result line output

**Release 2.00**

Released 5.8.97

General:

Re-designed for PFZ09 system

1024 x 768 half-screen layout to cooperate with Zymark software system

Windows NT 4.0

Networking features for data storage, parameter read/write and result output

Automatic:

On-line spectrum display to display background, standard and sample spectra in energy, abs or 2<sup>nd</sup> deriv.

Redesigned communication with Zymark system, result data written to common result file directly after being produced.

Removed auto-print

Manual:

Spectra set load/save included

Std abs, sample abs and % diss displayed

Two spectra graphs can be individually configured to display energy, abs or 2<sup>nd</sup> deriv.

**Release 1.54**

released 20.6.96

General:

-implemented/corrected ASPECT-equivalent 2nd deriv algorithm

-absorbance spectra Standard/Sample are now interpolated to equidistant Lambdas

-INI parameter wavelength range restriction now functional

Automatic:

Background Absorbance now checked as |Background|<?|MaxAllowedAbs|

Plot got new Axis Font

Manual:

Spectra display now shows absorbance Standard/Sample spectra

2nd deriv display now computed according to ASPECT

**Release 1.53**

released 13.5.96

Bug Fixed in 2-pages plot in case of 12 vessels (labels, data)

removed measuring times from headline of sample measurement table in report file

**Release 1.52**

released 28.2.96

New installation routine

**Release 1.51**

released 28.2.96

Std Abs in record file is now 6 dec places for deviation method

**Release 1.5**

released 1/96

rearranged Output &amp; Print sequence to avoid Application error (due to system overload)

**Release 1.4**

released 22.12.95

Manual Program:

Current PAR file now actualized in ZDIS\_MUX.INI

Automatic Program:

Added Error Message if PAR file is missing

General:

Included ASPECT plus type second deriv (Savitzky-Golay smooth, then SG9 deriv)

**Release 1.3**

released 20.12.95

Manual Program:

conc reading of sample to 5 decimal places

Quit button works ( DEBUG=0 in ZDIS\_MUX.INI !!)

Automatic Program:

Report plots (Diss. curves) separately for each bath

Report plot fixed value range is now 110% Diss

Mean %Diss in report is now down to 2 decimal places

Bug fixed in 2nd deriv computation on Use2ndDeriv

General:

Shutter closes after InitDevice

Fixed bug in ASPECT file storage (raw data)

Created LabView Smooth/Deriv function, avoiding the SMODER.DLL call. Now no DLL foreign to LabView is called.

New features:

Plot of spectra, derivs, calibration curve and measuring parameters

added Cmd Info to Status Bar

Report plot includes additional (0,0) starting point

Spectral range can be limited by means of new ZDIS\_MUX.INI parameters MaxLambda, MinLambda ([Device])

**Release 1.21**

released 06.12.95

fixed bug in checkbox labelling

set precision of std conc (manual) to 5

fixed bug in Read Float Numbers ( confusion with English/ German Windows decimal

***Release 1.2***

released 03.12.95

separator character), affected reading of floats in INI files  
added Cmd Info to Robot Simulation Prompt Boxes

changed dissolution graph checkbox array to 12  
added extended report according to Pfizer  
added Dissolution plot  
implemented new driver (MUX Flag etc.)  
added new entry VesselVol to ZY.INI  
changed Concentration reading to %Dissolution

***Release 1.1***

released 22.11.95

Created on-site at ZYMARK, 18.-22.11.95  
Complete rearrangement of software.  
Now two separate programs Manual & Automatic: DT\_EXE.EXE, MANU\_EXE.EXE  
Major changes in ZY.INI, ZDIS\_MUX.INI etc. ( see Manual R. 1.1 )

***Release 1.0***


released 10.10.95

Test release, did not meet specification due to information losses  
Sent to ZYMARK beginning of Nov 95

# Supplements

## Pinout of MCS 500 I/O ports

(MCS 500 Rear View)



Pin	In	Out
1	GND	GND
2	Dig IN 1	Dig OUT 1
3	Dig IN 2	Dig OUT 2
4	Dig IN 3	Dig OUT 3
5	Dig IN 4	Dig OUT 4
6	+ 5V	+ 5V

## DIO connections

Connector	Pin	on (low)	off (high)
DIO Out	2	Ready	Busy
	3	BI/Std OK	not OK
	4	Sample OK	not OK
	5	Hardware OK	not OK
	DIO In	2	Read new ZY.INI
3		Read BI/Std	
4		Timer start	
5		Read Sample	



## Algorithm for smoothing and derivatives according to Savitzky-Golay

Smoothed spectra and their derivatives are computed according to the Savitzky-Golay algorithm. This fast numerical procedure provides the new data points using a set of predefined coefficients.

$$y_{new_j} = \frac{1}{(\Delta X)^{order\_of\_deriv}} * \frac{1}{coeff_{n+1}} * \sum_{i=1}^n (y_{old_{j+i-(n+1)/2}} * coeff_i)$$

The spectra to be processed must be of equidistant. type, i.e. the x-data points must be evenly spaced by  $\Delta X$ . The standard SG coefficients are computed for a  $\Delta X$  of 1.

The coefficients have been computed according to R. Ludwig [see "Methoden der Fehler- und Ausgleichsrechnung"; VEB Deutscher Verlag der Wissenschaften, Berlin 1969 ].


Coefficients for smoothing,  
using  $n$  points

Coefficients for derivatives  
( $n=9$ )

i	n=5	n=7	n=9	n=11	n=13	n=17	n=21	n=25	i	1st	2nd	3rd	4th
1	-3	-2	-21	-36	-11	-21	-171	-253	1	-4	28	-14	14
2	12	3	14	9	0	-6	-76	-138	2	-3	7	7	-21
3	17	6	39	44	9	7	9	-33	3	-2	-8	13	-11
4	12	7	54	69	16	18	84	62	4	-1	-17	9	9
5	-3	6	59	84	21	27	149	147	5	0	-20	0	18
6	35	3	54	89	24	34	204	222	6	1	-17	-9	9
7		-2	39	84	25	39	249	287	7	2	-8	-13	-11
8		21	14	69	24	42	284	342	8	3	7	-7	-21
9			-21	44	21	43	309	387	9	4	28	14	14
10			231	9	16	42	324	422	10	60	462	198	143
11				-36	9	39	329	447					
12				429	0	34	324	462					
13					-11	27	309	467					
14					143	18	284	462					
15						7	249	447					
16						-6	204	422					
17						-21	149	387					
18						323	84	342					
19							9	287					
20							-76	222					
21							-171	147					
22							3059	62					
23								-33					
24								-138					
25								-253					
26								5175					

The described algorithm is used in the ZEISS software to produce smoothed spectra and their derivatives.

## Labview Graph Handles




The graph palette has controls for panning (scrolling the display area of a graph) and for zooming in and out of sections of the graph.


If you press the x **autoscale button**, LabVIEW autoscales the X data of the graph. If you press the y autoscale button, LabVIEW autoscales the Y data of the graph. If you want the graph to autoscale either of the scales continuously, click on the lock switch, shown at the left, to lock autoscaling on.




The **scale format buttons** give you run-time control over the format of the X and Y scale markers respectively.

You use the remaining three buttons to control the **operation mode** for the graph.

 Normally, you are in **standard** operate mode, indicated by the plus or crosshatch. In operate mode, you can click in the graph to move cursors around.

 If you press the **panning tool**, shown to the left, you switch to a mode in which you can scroll the visible data by clicking and dragging sections of the graph.

If you press the **zoom tool**, shown at the left, you can zoom in on a section of the graph by dragging a selection rectangle around that section. If you click on the zoom tool and hold for a moment, you get a pop-up menu you can use to choose some other methods of zooming. This menu is shown in the following illustration.



A description of each of these options follows.



Zoom by rectangle.



Zoom by rectangle, with zooming restricted to x data (the y scale remains unchanged).



Zoom by rectangle, with zooming restricted to y data (the z scale remains unchanged).

Undo  
Zoom

Undo last zoom. Resets the graph to its previous setting.



Zoom in about a point. If you hold down the mouse on a specific point, the graph continuously zooms in until you release the mouse button.



Zoom out about a point. If you hold down the mouse on a specific point, the graph continuously zooms out until you release the mouse button.

*Note: For the last two modes, zoom in and zoom out about a point, shift clicking zooms in the other direction.*

---

## ***Zeiss ASPECT plus data file format***

author: E.Stein  
released : 08.06.95

### ***General remarks***

- The header of the data file is similar to a windows-INI-file. To read a data file it is necessary to interpret the sections GENERAL and DATA.  
Every entry has to start at the beginning of a line. The section keywords are enclosed in square brackets ('[' and ']') and immediately followed by a CR/LF.
- The keywords in the sections must be followed immediately by an equal sign and then by the corresponding value of the respective keyword (in GENERAL section). In the DATA section the equal sign is to be followed immediately by a CR/LF. Here the first byte of the data has to stand on the beginning of the next line. After the last byte of data a CR/LF is to be appended.
- All data are 4-Byte IEEE float values. ASPECT will accept values in the range (-1.0 E15 .... 1.0 E15) only.
- The section GENERAL is always the first section in the data file.
- Several more sections and keywords are defined for internal use and may occur in the data files.
- The data file should be opened in binary mode.

---

## **Format**

### **[GENERAL]**

TITLE= small note of up to 50 characters

TYP= 0 old DOS-aspect format,  
1 old DOS-aspect format,  
2 old DOS-aspect format,  
3 old DOS-aspect format,  
4 current data file format

TIME= date and time as a string in the format:  
dd.mm.yyyy/hh:mm:ss

XUNITS= unit of the abscissa values

YUNITS= unit of the ordinate values

ZUNITS= unit of the cycle values

FIRSTX= start value of measurement range as string

LASTX= end value of measurement range as string

NPOINTS= number of measurement values per cycle as string

NCYCL= number of cycles as string

### **[DATA]**

ZDATA= cycle values in a field of float values (data type float). The size of the field (number of elements) is specified by keyword NCYCL. ZDATA is not necessarily included in the data file, even if the number of cycles is greater than 1 !

XDATA= abscissa values in a field of float values. The size of the field is specified by keyword NPOINTS. The abscissa values are the same for all cycles.

YDATA= measurement values in a field of float values. The size of the field is NPOINTS \* NCYCL. The cycle values are stored in subsequent order (cycle 1, cycle 2, ...)